





SMZ161 Series

SMZ-161-TP





SMZ161 Series

SMZ-161-TP

Optical system	Greenough
Observation tube	Trinocular head
Inclination	45° inclined
Trinocular light split	Fixed 50:50
Interpupillary distance	50-75mm
Diopter adjustment	On both tubes, +/- 5 diopter
Eyepieces	Widefield WF10X/20mm
Objectives system	Zoom, ratio 6:1
Objectives	0.75X~4.5X
Working distance	110mm
Stand type	Pole type stand without illumination
Head holder	For Ø25mm column and Ø76mm head
Focus mechanism	Coarse focusing system with tension adjustment
Focusing stroke	43mm
Accessories included	Dust cover
Dimensions LxWxH	280x220x367mm
Base	280x220mm
Column/Arm	Height 221mm, Ø25mm
Net weight	3.4kg